



**Tablets
Capsules
LCT's
Softgels**

Up to 60,000 PPH

Delta C02 Laser

One or Two-Sided C02 Laser Writer
Single-Lane Carrier Link Feed System

Delta C02 Laser

C02 Laser Marking for Tablets & Capsules

The Delta C02 Laser writer can laser write an alphanumeric, trademark logo, or symbol, soft gelatin capsule, or caplet. The fully adjustable hopper allows the product to flow smoothly to the laser station and can laser mark 50,000 to 70,000 pieces per hour.

Features

- Single or two-sided laser marking
- Adjustable hopper for maximum fill rate
- Servo drive system for precision print registration
- Quick change tooling
- Swingarm 15" touchscreen HMI
- Wi-Fi connectivity
- Batch recording
- 21 CFR Part 11 Compliant
- CE and UL certified

Capabilities


- Output Up to 60,000 / Hour
- Single or Double Sided Marking with independent graphics for each side
- Mark alpha-numeric text, vector graphics, logos and bar codes
- Graphics editing software provided
- Adjustable parameters to control logo size, line width and etch depth
- Mark Bi-Color products with logo oriented on selected color side
- Orient logos based on product feature, such as tablet bisect line

Options

- Vision inspection system with single piece rejection
- Validation document packages available
- Floor level infeed bulk conveyor

Safety & Regulatory Compliance

- | | |
|-----------------------------|--------------------|
| • cGMP for Pharmaceuticals: | • NFPA 79 |
| FDA 21 CFR Compliance: | • IEC / EN 60204-1 |
| - Part 210 | • ISO 13849-1 |
| - Part 211 | • CRC, c 870 |
| - Part 11 | • CDRH Compliant |
| • CE Compliant | |



Specifications



Made in the USA

H	2150mm / 84.6in
H1	1264.3mm / 49.8in
H2	1578.9mm / 62.2in
W1	829.2mm / 32.6in
W2	2416.4mm / 95.1in
L1	1663.5mm / 65.5in
L2	2328.5mm / 91.7in
POWER	240 VAC / 20 amp / 1 phase
AIR PSI	Air 90 psi @ 5 scfm